

**OXOID QUALITY ASSURANCE**

**PRODUCT SPECIFICATION**

**TRYPTONE SOYA AGAR**

**Typical Formula, Oxoid CM131**

**PO0193C**

	grams per litre
Tryptone	15.0
Soya peptone	5.0
Sodium chloride	5.0
Agar	15.0

**Preparation**

Suspend Tryptone Soya Agar (40g/litre) in deionised water. Sterilise at 121°C for 15 minutes. Aseptically dispense into contact plates. Label plates, wrap and label pack.

**Format**

Ten 55mm contact plates, wrapped in a single cellulose-based film wrap. Each plate is ink-jet printed with (abbreviated) product name, product code, lot number and expiry date.

**Labels**

Label gives details of product name, product code, recommended storage temperature, lot number and expiry date.

**Physical Characteristics**

Physical Tests

pH	- 7.3 ± 0.2
Colour	- Straw
Clarity	- Clear
Fill volume	- 17.75 ± 1.25ml

Packaging and presentation:

General appearance of pack and label should be satisfactory. Label data should be correct.

**Sterility Test**

Macroscopic examination should show no evidence of microbial growth after incubation at 20-25°C and 30-35°C for 5 days.

## Microbiological Testing Using Optimum Inoculum Dilution

### Results after incubation at 30-35°C for 2 days

#### Inoculum less than 100 colony forming units

<i>Staphylococcus aureus</i>	ATCC 6538	Straw colonies
<i>Escherichia coli</i>	ATCC 8739	Cream colonies
<i>Bacillus subtilis</i>	ATCC 6633	Irregular, straw colonies
<i>Pseudomonas aeruginosa</i>	ATCC 9027	Straw colonies

### Results after incubation at 20 - 25°C for 5 days

#### Inoculum less than 100 colony forming units.

<i>Candida albicans</i>	ATCC 10231	Cream colonies
<i>Aspergillus niger</i>	ATCC 16404	White mycelium, black spores

Colony counts shall be equal to or greater than 50% of the control medium.

#### Storage conditions

Store away from the light at between 8 – 15°C.