



MICROBACT™

IDENTIFICATION KITS



MICROBACT™ 2009
IDENTIFICATION PACKAGE V2.04
(WINDOWS™)

MB1244A

Oxoid Limited, Wade Road, Basingstoke, Hampshire RG24 8PW, UK

Phone: +44 (0) 1256 841144 Fax: +44 (0) 1256 463388

© Oxoid Australia Pty Limited, 2009



MICROBACT™ COMPUTER AIDED IDENTIFICATION SYSTEMS

FOR WINDOWS™ 95/98/2000/XP

CD INSTALLATION INSTRUCTIONS

(Please remove earlier versions of Microbact™ software before installation)

1. Exit all Windows programs
2. Insert the CD-Rom into the CD Drive
3. Click on the Start button. Choose Run from the File Menu.
In the Command Line Box type the letter of the CD drive followed by a colon, back slash, and the word setup and then press the Enter key.
eg. type D:\setup — press Enter
4. The installation program will be loaded
5. Answer all questions when prompted

Note: The Setup program can also be run by browsing the CD with Windows Explorer and double-clicking on the setup.exe file.