



Anaerogen W-Zip Compact AN0010W

For increased convenience, the new Anaerogen W-Zip Compact (AN0010W) includes everything you need for anaerobic atmosphere generation.

Oxoid Atmosphere Generation System (AGS) products quickly create the required anaerobic or micro-aerophilic conditions without generating hydrogen, safely facilitating prompt identification of cultures by maximising recovery and colony size. AnaeroGen™ reacts rapidly on contact with air to produce an atmosphere of <1% oxygen supplemented with carbon dioxide within 30 minutes, providing ideal conditions for growth of fastidious and obligate anaerobes.

The W-Zip pouches are transparent with an integral seal. The W-Zip is easy to use; simply pinch the seal together at one end and squeeze all the way across, ensuring there are no gaps.

No need for clips. Just add an indicator strip, before sealing.

AN0010W	AnaeroGen W-Zip Compact. Consta de 10 saquitos anaerogen compact y 10 bolsas de cierre en W-Zip.	10 sobres + 10 bolsas
----------------	--	-----------------------------

